



SURGICAL MASKS TECHNICAL DATA SHEET




- Fluid resistant;
- Comfortable, & hypoallergenic;
- Ultrasonic Sealed;
- +95% filtration;
- Latex free;
- Adaptable nose bar for better fitting;
- BFE Tested EN 1492001;
- 50 pieces per box.

Microbiologically cleaned with UV ray before packaging.




A GUIDE TO MASK USE




Before putting on a mask, clean hands with alcohol-based hand rub or soap and water.



Cover mouth and nose with mask and make sure there are no gaps between your face and the mask.



Replace the mask with a new one as soon as it is damp and do not re-use single-use masks.



To remove the mask: remove it from behind (do not touch the front of mask); discard immediately in a closed bin; clean hands with alcohol-based hand rub or soap and water.